



Plastalgin is ideal for:

- Case study models
- Orthodontic models
- Opposing models
- Provisionals
- Mouthguards
- Bleaching trays

Plastalgin — So quick, so easy, so clean, it makes a lasting impression.

Plastalgin is a popular multi-purpose, cost-effective alginate impression material that works so quickly and so easily that dental professionals are consistently impressed. And Plastalgin's mixed creamy homogenous consistency without excessive flow gives your patients a level of comfort not usually found with other alginate impression materials.



Materials

Some lead... others follow.

Making Impressions Quick and Easy...

Plastalgin is a highly versatile alginate impression material that is simple to use, easy to handle, and available in 2 formulations to deliver the perfect balance of working and setting times. Whether you choose Plastalgin Normal Set (Apple flavored) or Plastalgin Ortho-Quick Set (Apricot flavored), both products are naturally hydrophilic and deliver precise impressions with an excellent reproduction of marginal details.

Making impressions safe and dust-free...

Plastalgin's dust-free formulation is specifically designed to reduce the amount of airborne particles dispersed to exposed individuals. And because it is dust-free, cleanup is faster and easier with an overall environment that is healthier for you and your patients.

Dental professionals rely on Plastalgin for:

- Case study models
- Orthodontic models
- Fabrication of casts for partial denture frameworks, removable orthodontic or occlusal splint appliances
- Mouthguards
- Provisional restorations
- Bleaching trays

Physical properties:

	Plastalgin Normal Set	Plastalgin Ortho-Quick Set
Mixing Time:	45 seconds	30 seconds
Total Working Time:	2:15	1:15
Setting Time:	3:10	1:45
Viscosity:	38mm	36mm
Reproduction of Details:	20m	50m
Deformation When Compressed:	12.5%	12.5%

Ordering Information

ITEM NUMBER	PLASTALGIN TYPE
01-C1000	PLASTALGIN ORTHO QUICK SET, 17oz. — APRICOT FLAVOR
01-C1100	PLASTALGIN NORMAL SET, 17oz. — APPLE FLAVOR



Materials

Visit www.septodontusa.com for more information.

800-872-8305 | septodontusa.com

Some lead... others follow.